



PROMOTIONS | 2026
DU 1^{ER} SEPTEMBRE AU 31 DÉCEMBRE

Offres de fin d'année
AUTOMNE HIVER

**L'EXCELLENCE
AU SERVICE DES
CABINETS
DENTAIRES**



DEPPELER
Signature of excellence

OFFRE 5+1 OFFERT
SUR LES INSTRUMENTS EN DESTOCKAGE

CONTACTEZ-NOUS !



Matériel premium, design moderne, décliné en plusieurs couleurs pour s'adapter à votre environnement. Avec Prodentalalpes, offrez à votre cabinet des équipements alliant performance, confort et esthétique.

NOS ATOUTS

- ✓ **Qualité et fiabilité** – Des produits sélectionnés pour leur durabilité et leur efficacité
- ✓ **Confort et ergonomie** – Pensés pour améliorer l'expérience des praticiens et des patients.
- ✓ **Service expert et réactif** – Un accompagnement personnalisé pour répondre à vos besoins.

CONTACTEZ-NOUS

✉ contact@prodentalpes.fr

☎ 06 01 21 16 35



prodentalpes.fr

ZI DES DRAGIEZ, 30 RUE DE GRÉBELIN,
74800 LA ROCHE-SUR-FORON

INSTRUMENTATION SOINS

Une spatule universelle ! Faites tout avec.

D AESTHETIC RESTORATION KIT

Developped in collaboration with Dr. Von Sontagh

The Dr. Von Sontagh Restoration Kit contains 4 instruments (including 1 new instrument) to perform all esthetic restorations with precision and efficiency



OP43+ADEP BC

Kit Von Sontagh+ADEP

- OP80C + ADEP VC - Customized plugger with visual guides
- CPR05C + ADEP O - Customized sharp spatula with ruler
- CK10 + ADEP J - Modeling instrument
- OP43 + ADEP BC - Multi-purpose versatile spatula
- CleaNext B5 - Sterilisation box (p.63)

Instruments available for sale in kits only




SINGLE-END EXPLORERS

Thanks to the unique characteristics of our raw material, our explorers have an unmatched lifetime. Each working end is hand-finished to guarantee a perfectly sharp tip.



DOUBLE-END EXPLORERS




Sharp
Resilient and flexible




DH2 & DH3
For pocket calculus detection
DH2= angulation Gracey 11/12
DH3= angulation Gracey 13/14



CARVERS




HZ22
Holz



TH27
Hirt



SM20




SM18




SM19


HZ22, TH27, SM20, SM18
Angulation left, right



SM17




ZI54A
Mühlemann




VD8
Wiland

VD8
Approximal




0.9 mm
HZ8
Holz




0.7 mm
HZ9
Holz

HZ8, HZ9
Box carver to remove the excess

D SPATULAS ROUNDED AND SHARP



D SPATULAS-PLUGGERS



SPATULAS

Blade thickness 0.22 mm

	MB4C Besek	CT X*	
1.6 mm Approximal			
	OP551C	CT X*	
1.0 mm			
	OP55C	CT X*	
1.6 mm			
	ZI66C	CT X*	
1.5 mm			
	OP3DC	CT X*	
1.8 mm			
	OP26C	CT X*	
2.6 mm			
	OP22C	CT X*	
3.0 mm			
	ZI62C	CT X*	
2.5 mm			
	OP6C Heidemann	CT X*	
2.5 mm			
	OP54C Heidemann	CT X*	
2.0 mm			
	OP58C	CT X*	
2.6 mm			
	OP84 Heidemann		
2.4 mm			
	OP40 Heidemann		
2.3 mm			2.6 mm

MODELLING INSTRUMENTS



MF2C
Fuzzi

CT X*



PJ4



MB3C
Besek

CT X*



HZ3C

CT X*



CPRO3C

CT X*



MV1X+ADEP
Vasant



MO1
Monari



TH28C
Hirt



HZ25C

CT X*



CK28BC
Carver

CT*



OP37



28BR29



BALL PLUGGERS

	BR086	
ø 0.8 mm		ø 0.6 mm
	27SBR27SS	
ø 1.0 mm		ø 0.8 mm
	27BR27S	
ø 1.3 mm		ø 1.0 mm
	26BR27	
ø 1.7 mm		ø 1.3 mm
	25BR26	
ø 2.3 mm		ø 1.7 mm
	24BR25	
ø 2.7 mm		ø 2.3 mm
	08BR06	
ø 0.8 mm		ø 0.6 mm
	10BR08	
ø 1.0 mm		ø 0.8 mm
	MB6C Besek	
ø 2.3 mm		ø 0.6 mm
	127C	
ø 1.3 mm		ø 1.0 mm
	126C	
ø 1.7 mm		ø 1.3 mm
	125C	
ø 2.3 mm		ø 1.7 mm
	124C	
ø 2.7 mm		ø 2.3 mm
	123C	
ø 3.2 mm		ø 2.7 mm

PLUGGERS



MB1C
Besek

CT X*



MB2C
Besek

CT X*




OP80C

CT X*



OP46C

CT X*



OP46DC

CT X*



OP80PC

CT X*



OP36MPC



OP36P



OP36L



OP27



OP27D



PLUGGERS

1.3 x 0.6 mm



TH10
Hirt



ø 0.7 mm

L 3.0 mm



ZI51A



L 2.0 mm

L 2.2 mm



OP91



ø 1.1 mm

ø 0.9 mm



OP35



ø 0.6 mm

SERRATED PLUGGERS



ø 1.2 mm



OP36M



ø 0.9 mm

ø 1.6 mm



OP36D



ø 1.1 mm

ø 2.2 mm



OP36



ø 1.6 mm

ø 1.7 mm



OP52



ø 1.0 mm

ø 2.6 mm



ZI51



ø 2.2 mm

ø 2.2 mm



OP44



ø 1.8 mm

ø 2.6 mm



1BR2



ø 1.7 mm

ø 3.6 mm



2BR3



ø 2.6 mm

M23®


Our mastery of the raw material not only gives the instruments their incomparable cut, but also an extremely hard-wearing.



«Working grip and ergonomic design are impressive
-Required force is spread evenly on the hand
-Steel quality is excellent, sharpness remains longer
-Tolerates well sharpening – just a featherlike « kiss » is sufficient

-Steel quality is extremely important for clinicians – and only he/she can feel it
-Deppeler angulation requires less wrist stressing movements
-High quality instruments play vital role for hand instrumentation»

Päivi Maunola, Helsinki, Finland




POSTERIOR


Premolars, Molars for approximal surfaces



Fine versions
These fine scalers enable extremely delicate debridement.
IMPORTANT: not only the blade, but the entire shaft is thinner. This allows greater sensitivity during the clinical work.
It is also a little less strong.




M23A®



Fine versions

These fine curettes enable extremely delicate debridement.
 IMPORTANT: not only the blade, but the entire shaft is thinner. This allows greater sensitivity during the clinical work.
 It is also a little less strong.


PREMOLARS, MOLARS




INSTRUMENTATION IMPLANTS & SINUS

IMPLANT POCKET PROBE IN PEEK

1. Adapts to the anatomy due to its flexibility
2. Indelible precise visual guide through the material
3. Also suitable for probing of natural teeth



DOUBLE-END PROBES




ENGRAVED PROBES



BIOTYPE PROBE

Biotype probe clinical use:

1. If the fine / graduated part is visible through the gingiva, then it is a very thin biotype
2. If the fine / graduated part is not visible but the thick / black part is, then it is a medium biotype
3. If the thick / black part is not visible through the gingiva, then it is a very thick biotype



DBS12 + ADEP B

Fischer, Kai R et al. "Gingival biotype revisited-novel classification and assessment tool." *Clinical oral investigations* vol. 22,1 (2018): 443-448.

Black side $\varnothing$ 0.7 mm
Graduated side (3-6-9-12) $\varnothing$ 0.45 mm



Dr. Kai Fischer, specialist dentist for periodontology (ZÄKWL), specialist for periodontology (DGParo e.V.) «Dental specialists at the cathedral».

University of Zurich,
University of Witten/
Herdecke

SINGLE-END PROBES

Deppeler periodontal probes combine flexibility and resilience. The pocket visual guides are very long-lasting. Each tip is perfectly rounded to ensure atraumatic probing.



Photo: L. Parisi, Italy



HH12 DMS-N



HH12
3-6-9-12
 $\varnothing$ 0.48 mm



HH11
3-6-8-11
 $\varnothing$ 0.48 mm



DHZ12
3-6-9-12
 $\varnothing$ 0.48 mm



UNC15
1-15
 $\varnothing$ 0.48 mm



GY12 Günay
1-2-3.5-5.5-8.5-11.5
 $\varnothing$ 0.50 mm



DHB11
3-6-8-11
 $\varnothing$ 0.50 mm



KHR12
3-6-9-12
 $\varnothing$ 0.50 mm



WHO
3.5-5.5-8.5-11.5
 $\varnothing$ 0.50 mm


EFFICIENT WITHOUT DAMAGING

These solid titanium curettes are softer than the titanium used for the production of implants.
Surface of a titanium implant under the microscope.*


new



treated with a steel curette



treated with a titanium curette



* L. Parisi, C. Lorenzi, A. Genovesi (2015)
"Trattamento di colli implantari con tre tipi di curette :
TITANIO - ACCIAIO - PLASTICA".


TITANIUM SMART SCALING

2 instruments for all clinical cases.



«The new designed Titanium instruments from Deppeler are a class of their own. Better access to tight soft tissue niches or areas around the superstructure is achieved due to the development of a smaller yet durable working end.»

Bernita Bush Gissler, RDH,BS, Bern, Switzerland



Titanium Smart Scaling kit | KITTISCN

TIS1CN
TIS2CN
AD23CN
PP12DMS-N
MIR-SDMS-N
CLEANext B5 - (p. 63)



Additional study :
A. Butera, M. Delgrosso, A. Genovesi, L. Parisi (2015) 'Efficacy of titanium instrumentation combined with air-polishing with glycine for the reduction of mucositis and peri-implantitis: clinical studies compared' Hygiene Tribune Italian edition 12:24-25.


NABER'S PROBES



FURCATION CURETTES

2 steps for clean furcations:

1. Curettage with the curettes F10 or F11 depending on the size of the furcation (see pictures).



2. Smoothing with diamond coated curettes F12 Nabers like diamond coated curette.
Very fine diamond grain not to damage the root surface.

> clinical study p.18




TUNNELING KNIVES

«The Deppeler tunneling knives are exceptionally well designed instruments, perfectly suited for minimalinvasive surgical interventions. The surgical set is ideal for tunnelised flap techniques aiming at root coverage, enlargement of keratinised tissues and socket preservation.»


Prof. Péter Windisch /Department of Periodontology, Semmelweis University, Budapest, Hungary

D.M.D., Semmelweis University Faculty of Dentistry, Budapest, Hungary, 1982-1987
Associate Professor /Department of Periodontology, Semmelweis University, Budapest, Hungary, 2005-2012
Professor /Department of Periodontology, Semmelweis University, Budapest, Hungary, 2012-present
Chairman, Hungarian Society of Periodontology
Vice chairman, Hungarian Society for Dental Implantology



Tunneling knives kit | KITTK+ADEP

TKL+ADEP BC - Tunneling knife left
TKR+ADEP BC - Tunneling knife right
TKF+ADEP BC - Tunneling knife papillae
TKI+ADEP BC - Tunneling knife anteriors
TKP+ADEP RS - Tunneling knife papillae
CLEANext B5 - Sterilization box - (p. 63)




More on www.deppeler.ch

SINUSLIFT



KIT SINUS LIFT | KITSINUS

CU3+ADEP R
AF30+ADEP O
CU5+ADEP J
SL2+ADEP V
BS1+ADEP BC
2BR3+ADEP VI
CLEANNext B10 - (p. 63)



Instruction for use of the SL360 curette

"After having reached the sinus membrane through osteotomy, the release or detachment of the membrane at the edge of the window will be facilitated due to the atraumatic elephant foot curette. Thanks to its geometric form, all sides and edges of the sinus window can be instrumented. In addition, removal of bone spicules will be possible."


Dr. med. dent. Cionca Norbert,
specialist in periodontology,
senior lecturer at the University of Geneva, Switzerland



PERIOTOMES

Fine and rounded blade.

The hand-finished cutting edge enables you to work efficiently and gently. These instruments guarantee you the best possible preparation for an extraction.



SYNDESMOTOMES

These instruments have been designed to cut ligaments before an extraction.

Their top quality manufacture and their exceptional cutting edge guarantee precise and efficient work.

IMPORTANT: these instruments should never be used as elevators.



The Deppeler straight-line desmotome VT1 is one of the best instruments for removing root tip debris. Throughout my professional career, I have tried new instruments to deal with this problem. Again and again I came to the conclusion that the Deppeler syndesmotom is the best possible instrument for this task. That is why I still use this instrument for this indication and can recommend it to any dentist practising surgery. I have also tested many modifications: bigger, smaller, thicker and also thinner. But I have to say it is perfectly balanced.

Prof. Dr. Markus Hürzeler, Hürzeler/Zuhr Dental Practice, Munich, Germany




SMART SCALING CURETTES

All root surfaces can be treated with only two instruments. Deppeler Smart Scaling's unique geometry, revolutionary angulation and double-edge blades, enables perfect access to periodontal pockets.



The Smart Scaling curettes are also available with 3-6-9-12 visual guides for maximum precision.

Anteriors & Buccal/Palatal Smart Scaling curette with 3-6-9-12 visual guides | DSS1CN
Mesial/Distal Smart Scaling curette with 3-6-9-12 visual guides | DSS2CN



WIDER ANGLE = SAFE TREATMENT

The Smart Scaling blade has a double edge with a wider angle and a rounded conical shape. This innovative geometry enables a safe debridement without damaging the healthy tissue.



Smart Scaling kit | KITDS(S)CN

- In this kit:
- D(S)S1CN
 - D(S)S2CN
 - M23CN®
 - H125
 - MIR-S DMS-N
 - CLEANNext B5 - (p. 63)



How to use the
Smart Scaling
curettes



"Modern hand instruments should meet the requirements of minimal invasive treatment. They need to be precise and delicate, whilst retaining strength. Besides ergonomic considerations, they should enable the most efficient clinical work. The double sided cutting edges of the Smart Scaling curettes and the special angulation of the working end make it possible to access all surfaces with only two instruments. The periodontal pocket visual guides on the end of the instruments enable a local precise debridement without damaging healthy tissue. The CLEANNext handle reduces operative fatigue. Try the Smart Scaling curettes, you will be enthusiastic!"

Daniela Krejci, Dental hygienist HF, CDS Signy, Switzerland





depuis 2011

L'EXCELLENCE AU SERVICE DES CABINETS DENTAIRES

CONTACTEZ-NOUS OU COMMANDEZ SUR
NOTRE SITE INTERNET POUR BÉNÉFICIER
DE CES OFFRES EXCLUSIVES.

PRODENTALPES

 06 01 21 16 35

 contact@prodentalpes.fr

 ZI DES DRAGIEZ, 30 RUE DE GRÉBELIN,
74800 LA ROCHE-SUR-FORON



10 000 PRODUITS TOUTES MARQUES

Totalement indépendant des laboratoires, nous sommes à votre écoute
et nous vous conseillons de façon totalement objective.

Visitez le site prodentalpes.fr



Découvrez nos Promotions & Formations